Hardware User's Manual

Temperature Control Unit

Homeothermic Blanket

References:

HB101 (76-0385), RS (76-0111)

Version:

1.0

Limitation of Liability

PANLAB does not accept responsibility, under any circumstances, for any harm or damage caused directly or indirectly by the incorrect interpretation of what is expressed in the pages of this manual.

Some symbols may have more than one interpretation by professionals unaccustomed to their usage.

PANLAB reserves the right to modify, in part or in total, the contents of this document without notice.

1. SYMBOLS TABLE

Recognising the symbols used in the manual will help to understand their meaning:

DESCRIPTION	SYMBOL
Warning about operations that must not be done because they can damage the equipment	
Warning about operations that must be done, otherwise the user can be exposed to a hazard.	
Protection terminal ground connection.	
Warning about a hot surface which temperature may exceed 65°C	
Warning about a metal surface that can supply electrical shock when it's touched.	
Decontamination of equipments prior to disposal at the end of their operative life	
Waste Electrical and Electronic Equipment Directive (WEEE)	

2. GOOD LABORATORY PRACTICE

Check all units periodically and after periods of storage to ensure they are still fit for purpose. Investigate all failures which may indicate a need for service or repair.

Good laboratory practice recommends that the unit be periodically serviced to ensure the unit is suitable for purpose. You must follow preventive maintenance instructions. In case equipment has to be serviced you can arrange this through your distributor. Prior to Inspection, Servicing, Repair or Return of Laboratory Equipment the unit must be cleaned and decontaminated.

Decontamination prior to equipment disposal

In use this product may have been in contact with bio hazardous materials and might therefore carry infectious material. Before disposal the unit and accessories should all be thoroughly decontaminated according to your local environmental safety laws.

3. UNPACKING AND EQUIPMENT INSTALATION

WARNING: Failure to follow the instructions in this section may cause equipment faults or injury to the user.

- A. No special equipment is required for lifting but you should consult your local regulations for safe handling and lifting of the equipment.
- B. Inspect the instrument for any signs of damage caused during transit. If any damage is discovered, do not use the instrument and report the problem to your supplier.
- C. Ensure all transport locks are removed before use. The original packing has been especially designed to protect the instrument during transportation. It is therefore recommended to keep the original carton with its foam parts and accessories box for re-use in case of future shipments. Warranty claims are void if improper packing results in damage during transport.
- D. Place the equipment on a flat surface and leave at least 10 cm of free space between the rear panel of the device and the wall. Never place the equipment in zones with vibration or direct sunlight.
- E. Once the equipment is installed in the final place, the main power switch must be easily accessible.
- F. Only use power cords that have been supplied with the equipment. In case that you have to replace them, the spare ones must have the same specs that the original ones.

- G. Make sure that the AC voltage in the electrical network is the same as the voltage selected in the equipment. Never connect the equipment to a power outlet with voltage outside these limits.

WARNING

For electrical safety reasons you only can connect equipment to power outlets provided with earth connections . This equipment can be used in installations with category II over-voltage according to the General Safety Rules.

The manufacturer accepts no responsibility for improper use of the equipment or the consequences of use other than that for which it has been designed.

PC Control

Some of these instruments are designed to be controlled from a PC. To preserve the integrity of the equipment it is essential that the attached PC itself conforms to basic safety and EMC standards and is set up in accordance with the manufacturers' instructions. If in doubt consult the information that came with your PC. In common with all computer operation the following safety precautions are advised.

- WARNING**
- To reduce the chance of eye strain, set up the PC display with the correct viewing position, free from glare and with appropriate brightness and contrast settings
 - To reduce the chance of physical strain, set up the PC display, keyboard and mouse with correct ergonomic positioning, according to your local safety guidelines.

4. MAINTENANCE

WARNING: Failure to follow the instructions in this section may cause equipment fault.

- **PRESS KEYS SOFTLY** – Lightly pressing the keys is sufficient to activate them.
- Equipments do not require being disinfected, but cleaned for removing urine, faeces and odour. To do so, we recommend using a wet cloth or paper with soap (which has no strong odour). **NEVER USE ABRASIVE PRODUCTS OR DISSOLVENTS.**
- **NEVER** pour water or liquids on the equipment.
- Once you have finished using the equipment turn it off with the main switch. Clean and check the equipment so that it is in optimal condition for its next use.
- The user is only authorised to replace fuses with the specified type when necessary.


Figure 1. Power inlet, main switch and fuse holder.

FUSE REPLACEMENT OR VOLTAGE SETTING CHANGE

In case of an over-voltage or other incident in the AC net making it impossible to turn on the equipment, or if the equipment voltage setting is incorrect, check fuses according to the following procedure.

- 1 Remove power cord from the power inlet.

- 2 Open fuse-holder by pulling the flange with a regular screwdriver.


Figure 2. Open fuse-holder door.

- 3 Extract fuse holder using the screwdriver.



Figure 3. Extract fuse-holder.

- 4 Replace fuses if necessary. Insert fuses in the fuse-holder in the correct position.

CORRECT

INCORRECT

Figure 4. Fuses position.

- 5 Insert the fuse-holder again, positioning it according to the voltage in the AC net.

115V POSITON

230V POSITION

Figure 5 Fuse holder position.

- 6 If the fuses blow again, unplug the equipment and contact technical service.

WARNING

For electrical safety reasons, never open the equipment. The power supply has dangerous voltage levels.

5. TABLE OF CONTENTS

1. SYMBOLS TABLE	2
2. GOOD LABORATORY PRACTICE	2
3. UNPACKING AND EQUIPMENT INSTALATION	3
4. MAINTENANCE	5
5. TABLE OF CONTENTS	7
6. INTRODUCTION	9
7. EQUIPMENT DESCRIPTION	10
7.1. FRONT PANEL	10
7.2. REAR PANEL	11
7.3. HOMEOTHERMIC BLANKET	12
8. EQUIPMENT CONNECTION	13
8.1. HB101/2	13
8.2. HB101/2 RS	14
8.3. CONNECTING THE BLANKET	15
8.4. CONNECTING THE TEMPERATURE PROBE	15
9. WORKING WITH THE EQUIPMENT	16
9.1. CONDUCTING AN EXPERIMENT	16
9.2. CLEANING THE TEMPERATURE PROBE	17
9.3. TEMPERATURE PROBE DISINFECTION	17
9.4. TEMPERATURE PROBE STERILIZATION	17
9.5. PAD OR BLANKET CLEANING	17
10. USING THE SEDACOM SOFTWARE (RS VERSION)	18
11. TROUBLESHOOTING	19

12.	PREVENTIVE MAINTENANCE	20
13.	SPECIFICATIONS	21

6. INTRODUCTION

The HB 101/2 is a device used for monitoring the warming and temperature of homeothermic blankets for mice, rats and rabbits. The goal is the maintenance of body temperature in experimental animals. A rectal probe monitors the animal's body temperature, and this value is shown on the **TEMPERATURE** display.

Figure 6. HB101/2.

The RS versions enhance the standard model with the capability to transmit the current probe temperature to a PC and save it for future reference, using the **Sedacom** program (to be purchased separately).

7. EQUIPMENT DESCRIPTION

7.1. FRONT PANEL


Figure 7. HB101/2 Front Panel.

- **SELECT:** This button is used to set the blanket temperature between 32.0°C and 45.0°C. The display will show a flashing decimal dot, and the temperature setting can be modified with the ARROWS. The display will show the blanket temperature again if SELECT is pressed, or if 15 seconds pass without any button being pressed.
- **DISPLAY:** The display shows the blanket temperature with a resolution of 0.1°C. The display will show "---" if the temperature probe is not connected or the temperature is lower than 20°C. The decimal dot flashes when the device is in temperature selection mode.
- **HEAT:** This LED remains on while the blanket is heating.
- **ARROWS:** The selected temperature can be modified with the arrows after pressing the SELECT button. The temperature range is between 32.0°C and 45.0°C.

7.2. REAR PANEL


Figure 8. HB 101/2 Rear Panel.

- **BLANKET:** This connector is used to plug in the homeothermic blanket.
- **PROBE:** This connector is used to plug in the temperature probe.
- **TEST:** This button is used to check the control unit. Unplug the probe and the blanket, then select a temperature of 40.5°C by pressing the SELECT button, followed by the ARROWS. Press SELECT again to confirm the selection, then press the TEST button. The display will show approximately 40.0°C and the "HEAT" LED will be on because the selected temperature is higher than the displayed temperature. Now select 39.5°C and press the TEST button. The display will show approximately 40.0°C and the "HEAT" LED will go off, because the selected temperature is lower than the displayed temperature.
- **POWER INLET:** Power inlet, main switch and fuse holder.
- **RS-232:** This is a 9 pin female delta connector only available in the RS version. It is used to send data to a PC for processing with the **Sedacom** software.

7.3. HOMEOTHERMIC BLANKET

Connect it to the Blanket Connector on the Rear Panel.

Figure 9. Homeothermic blanket.

The homeothermic blanket models are as follows:


DESCRIPTION	DIMENSIONS (mm)
Soft blanket for rabbit	400x700
Hard blanket for rat	105x230
Soft blanket for rat	150x250
Hard blanket for mouse	90x114
Soft blanket for mouse	100x150

8. EQUIPMENT CONNECTION

8.1. HB101/2

The necessary connections are shown in the following schematic.

Figure

10.

Equipment connections.

The cables and connections are listed in the following table:

	FROM	TO	CABLE
1	HB101/2 blanket	Blanket	Blanket cable
2	HB101/2 probe	Animal rect	Temperature Probe

8.2. HB101/2 RS

HB101/2 RS has an RS-232 port to send data to the computer for processing with the **Sedacom** software. The necessary connections are shown in the following schematic.

Figure 11. Equipment connections.

The connections and cables are listed in the following table.

	FROM	TO	CABLE
1	HB101/2 blanket	Blanket	Blanket cable
2	HB101/2 probe	Animal rect	Temperature Probe
3	HB101/2 RS232	PC Com port	RS232 cable

8.3. CONNECTING THE BLANKET

To connect the blanket, align the connector with the socket (the mark must be visible on the upper part of the connector), insert the connector and then turn the metal ring to fasten it.


Figure 12. Connecting the blanket to the control unit.

8.4. CONNECTING THE TEMPERATURE PROBE

Insert the jack into the rear panel to connect the temperature probe to the control unit.



Figure 13. Connecting the probe to the control unit.

9. WORKING WITH THE EQUIPMENT


9.1. CONDUCTING AN EXPERIMENT

- 1 Connect the blanket and the temperature probe to the control unit rear panel.
- 2 Connect the power cord and turn the control unit on.


- 3 The display will show the probe temperature. If the probe temperature is lower than 20°C, or the temperature probe is not connected, or the temperature probe is damaged, the display will show "---".

- 4 Press the SELECT button to select the temperature. The decimal dot will flash. Select the temperature with the ARROWS. Press SELECT again to save the changes, or wait 15 seconds without pressing any button. The display will then show the probe temperature and the decimal dot will stop flashing.

- 5 Prepare the animal in the blanket. Place the probe in its anus.
- 6 Now, the control unit will work to reach the target temperature. The HEAT led will stay on while the blanket is heating.

- 7 The TEST button is to ensure that the unit is working properly. Disconnect the probe and push this button. A value of approximately 40.0°C should be displayed. If the selected temperature is lower than 40.0°C the HEAT led will not come on. If the selected temperature is higher than 40.0°C the HEAT led will come on.

The temperature probe must always be in contact with the animal or the blanket. If the probe is in the ambient air, the blanket may exceed the pre-set temperature, potentially over-warming or injuring the animal.

9.2. CLEANING THE TEMPERATURE PROBE

Read the temperature probe manual in order to follow the manufacturer recommendations.

9.3. TEMPERATURE PROBE DISINFECTION

Read the temperature probe manual in order to follow the manufacturer recommendations.

9.4. TEMPERATURE PROBE STERILIZATION

Read the temperature probe manual in order to follow the manufacturer recommendations.

9.5. PAD OR BLANKET CLEANING

To clean the pad or the blanket you can use a lightly wet cloth and then dry them with a dry cloth. If they're too dirty you can wet the cloth with a soapy solution to clean them, then remove foam with a wet cloth and finally dry them with a dry cloth.

10. USING THE SEDACOM SOFTWARE (RS version)

The purchase of the **Sedacom** software option is needed for transferring the data to a computer (please contact your local provider for more information). The **Sedacom** software reference is composed a USB Flash key containing the software Installer, License for use and **Sedacom** User's Manual. Follow the next instructions:

- Please refer to the **Sedacom** User's Manual for instructions on how to install and use the software with the present device.
- A serial port (RS232) communication cable (provided with the present device) is needed for the connection of the present device to the computer in which the **Sedacom** software is installed. Please refer to the Connection chapter 8 of the present User's Manual for instructions on how to connect this cable to the device.
- If the computer doesn't have any serial port, the RS232/USB adapter is needed (ref. CONRS232USB, contact your local provider for more information)

11. TROUBLESHOOTING

This table features instructions to solve most frequent problems.


PROBLEM	SOLUTION
The equipment does not start up.	<ul style="list-style-type: none"> • Ensure that the voltage of mains is the same as that selected in the fuse holder. • Check the condition of the fuses.
The display of the HB101/2 shows "---" instead of a numeric value.	<ul style="list-style-type: none"> • Ensure that the temperature probe is connected to the control unit. • With the help of an auxiliary thermometer check that the blanket temperature is not lower than 20°C or higher than 47°C (measuring range of the control unit). • Press the TEST button on the rear panel, if the display shows a temperature around 40.0°C and when you release the button the display shows "---", this means that the temperature probe is damaged and must be replaced by a new one.
The blanket does not heat	<ul style="list-style-type: none"> • Check that the blanket is connected to the control unit. • Check that the temperature probe is connected to the control unit. • Check that selected temperature is higher than the current temperature.
The blanket overheats above the selected temperature.	<ul style="list-style-type: none"> • Some overshoot of the temperature is expected. Wait a while to stabilize, if it returns to the selected temperature, the machine is working properly. If warming continues indefinitely despite the led HEAT being off contact the technical service.
The blanket feels hot but the temperature on the display varies slightly and does not reach the selected value.	<ul style="list-style-type: none"> • Ensure that the temperature probe remains in contact with the blanket or rectum of the animal whose temperature is regulated.

12. PREVENTIVE MAINTENANCE

	EXPERIMENT	MONTHLY
PAD OR BLANKET CLEANING	<input checked="" type="checkbox"/>	
CLEAN TEMPERATURE PROBE	<input checked="" type="checkbox"/>	
DISINFECT TEMPERATURE PROBE	<input checked="" type="checkbox"/>	
STERYLIZE TEMPERATURE PROBE		<input checked="" type="checkbox"/>
CHECK TEMPERATURE PROBE PLACING	<input checked="" type="checkbox"/>	

13. SPECIFICATIONS

POWER SUPPLY Input voltage: Frequency: Fuse: Maximum power: Conducted noise: Equipment start-up time:	115 /230 VAC 50 /60 Hz 2 fuses 5x20mm 1A 250V Fast 60 W EN55022 /CISPR22/CISPR16 class B <1 min
WARMING SPECIFICATIONS Temperature range: Display resolution: Accuracy: Highest temperature: Maximum output voltage in blanket: Minimum voltage in blanket: Noise: Noise frequency:	32–45 °C 0.1 °C +/-0.3 °C 65 °C (without probe feedback) 140 VDC 30 VDC 200 mV (peak to peak) 200 Hz
TEMPERATURE PROBE Technology: Measurement range showed: Linearity: Accuracy:	Thermistor 20-47 °C +/-0.1 °C +/- 0.1 °C
ENVIRONMENTAL CONDITIONS Operating temperature: Operating relative humidity: Storage temperature:	10°C to +40°C 0% to 85% RH, non-condensing 0°C to +50°C, non-condensing
COMUNICATIONS OUTPUT Standard interface: Connector:	RS232C Delta 9 contacts female connector
DIMENSIONS Width x Height x Depth: Weight:	161 x 115 x 250 mm 2.5 kg

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(GB) Note on environmental protection:

After the implementation of the European Directive 2002/96/EU in the national legal system, the following applies:

Electrical and electronic devices may not be disposed of with domestic waste. Consumers are obliged by law to return electrical and electronic devices at the end of their service lives to the public collecting points set up for this purpose or point of sale. Details to this are defined by the national law of the respective country. This symbol on the product, the instruction manual or the package indicates that a product is subject to these regulations. By recycling, reusing the materials or other forms of utilising old devices, you are making an important contribution to protecting our environment.

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